Excelsior[®] XT-17[™]

Microcatheter

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	10
Gebbruiksaanwijzing	12
Instruções de Utilização	14

Excelsior[®] XT-17[™]

Microcatheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Stryker Neurovascular representative.

For single use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

Stryker Neurovascular Excelsior XT-17 Microcatheters are single lumen devices designed to aid the physician in accessing distal vasculature when used with a guide catheter and steerable guidewire. Graded shaft stiffness ranging from a highly flexible tip to a semi-rigid proximal section aids the physician in tracking over selectively placed guidewires. A luer fitting located on the microcatheter hub is used for the attachment of accessories. A radiopaque tip facilitates fluoroscopic visualization.

Stryker Neurovascular hydrophilically coated Excelsior XT-17 Microcatheters are coated on the outer surface with Hydrolene[®] Coating that reduces friction during manipulation in the vessel.

Contents:

- One (1) Microcatheter
- One (1) Packaging Mandrel (for Preshaped Microcatheters only)
- One (1) Steam Shaping Mandrel
- One (1) Peel Away Introducer
- One (1) Steam Shaping Instruction Card

Table 1. Compatibility Information

Refer to the corresponding Stryker Neurovascular occlusion device Directions for Use for information regarding microcatheter compatibility.

Microcatheter	Effective Length (cm)	Distal Shaft Length (cm)	Microcatheter ID mm (in)	Min Guidewire ID mm (in)	Max Guidewire Dia. mm (in)	Non-Fibered Occlusion Coil Dia. mm (in)
Excelsior XT-17	150	7.5	0.43 (0.017)	1.0 (0.038)	0.36 (0.014)	0.24-0.38 mm (0.0095 in- 0.015 in)
Excelsior XT-17 Flex	150	15	0.43 (0.017)	1.0 (0.038)	0.36 (0.014)	0.24-0.38 mm (0.0095 in- 0.015 in)

INTENDED USE/INDICATIONS FOR USE

Stryker Neurovascular's Excelsior XT-17 Microcatheters are intended to assist in the delivery of diagnostic agents, such as contrast media, and therapeutic agents, such as occlusion coils, into the peripheral, coronary and neuro vasculature.

CONTRAINDICATIONS

None known.

WARNINGS

These devices are intended for use only by physicians trained in performing endovascular procedures.

- The accessories are not intended for use inside the human body.
- Limited testing has been performed with solutions such as contrast media, saline and suspended embolic particles. The use of these microcatheters for delivery of solutions other than the types that have been tested for compatibility is not recommended. Do not use with glue or glue mixtures.
- Carefully inspect all devices prior to use. Verify shape, size and condition are suitable for the specific procedure.
- Exchange microcatheters frequently during lengthy procedures that require extensive guidewire manipulation or multiple guidewire exchanges.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the microcatheter or guidewire against resistance could dislodge a clot, perforate a vessel wall, or damage microcatheter and guidewire. In severe cases, tip separation of the microcatheter or guidewire may occur.

PRECAUTIONS

- To reduce the probability of coating damage in tortuous vasculature, use a guide catheter with a minimum internal diameter as specified in Table 1 above, and is recommended for use with Stryker Neurovascular hydrophilically coated microcatheters.

- To control the proper introduction, movement, positioning and removal of the microcatheter within the vascular system, users should employ standard clinical angiographic and fluoroscopic practices and techniques throughout the interventional procedure.
- Exercise care in handling of the microcatheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- Use the product prior to the "Use By" date printed on the label.
- Limited testing indicates that Excelsior XT-17 Microcatheter is compatible with Dimethyl Sulfoxide (DMSO). The compatibility of Excelsior XT-17 Microcatheter with individual agents suspended in DMSO has not been established.

ADVERSE EVENT INFORMATION

Potential adverse events associated with the use of microcatheters or with the endovascular procedures include, but are not limited to:

- Access site complications
- Allergic reaction
- Aneurysm perforation
- Aneurysm rupture
- Death
- Embolism (air, foreign body, plaque, thrombus)
- Hematoma
- Hemorrhage
- Infection
- Ischemia
- Neurological deficits
- Pseudoaneurysm
- Stroke
- Transient Ischemic Attack
- Vasospasm
- Vessel dissection
- Vessel occlusion
- Vessel perforation
- Vessel rupture
- Vessel thrombosis

Adverse Event Reporting

Please notify your Stryker Neurovascular representative immediately if a device malfunctions or patient complication or injury is experienced or suspected. Please make every attempt to retain any suspect device, its associated components and their packaging for return to Stryker Neurovascular.

HOW SUPPLIED

Stryker Neurovascular products are sterile and non-pyrogenic in unopened packaging that is designed to maintain sterility unless the primary product pouch has been opened or damaged.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Required Additional Items

- Two 3 cm³ (3 cc) syringes and saline flush bags.
- Continuous flush setup with two rotating hemostatic valves (RHVs) (Tuohy-Borst type).
- Two stopcocks.

Preparations For Use

Warning: Inspect product before use for any bends, kinks or damage. Do not use a microcatheter that has been damaged. Damaged microcatheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

Caution: Wet dispenser coil or packaging tray and hydrophilically coated outer shaft of microcatheters prior to removal from packaging tray. Once the microcatheter has been wetted, do not allow to dry.

1. Gently remove the microcatheter from the dispenser coil or packaging tray and inspect the microcatheter prior to use to verify that it is undamaged.
2. If using a pre-shaped microcatheter, remove packaging mandrel from microcatheter tip and discard mandrel.

Caution: The packaging mandrel is not intended for reuse. The packaging mandrel is not intended for use inside the human body.

3. If desired, a shape proximal to the distal tip may be formed by following the steam shaping instructions provided.
4. Connect flush-filled syringe to microcatheter hub and flush microcatheter lumen to facilitate guidewire insertion.

Continuous Flush

The recommended continuous flush set up is illustrated in Figure 1. The RHVs provide a fluid tight seal and are attached to the guide catheter and microcatheter. Attach the stopcocks to the RHV sidearms, which become infusion ports for appropriate flush or contrast media injection.

Caution: Check that all fittings are secure so that air is not introduced into guide catheter or microcatheter during continuous flush.

Caution: In order to achieve optimal performance of Stryker Neurovascular Microcatheters and to maintain the lubricity of the Hydrolene® Coating surface, it is critical that a continuous flow of appropriate flush solution be maintained between the Stryker Neurovascular Microcatheter and guide catheter, and the microcatheter and any intraluminal device. In addition, flushing aids in preventing contrast crystal formation and/or clotting on both the intraluminal device and inside the guide catheter and/or the microcatheter lumen.


Figure 1. Example of Continuous Flush Setup

- Connect both RHVs to flush solution pressure bag. Pressure packs are convenient for this purpose. Maintain pressure above arterial pressure (approximately 40 kPa (300 mm/Hg)) in order to prevent retrograde flow of blood into guide catheter or microcatheter.
- Set up continuous flush by connecting first RHV to guide catheter and second RHV to hub of microcatheter.
- Carefully remove guidewire from its packaging and prepare the guidewire according to manufacturer's instructions.
- Carefully insert guidewire into microcatheter hub either directly or via guidewire introducer. Advance through microcatheter lumen until guidewire extends from distal microcatheter tip. If guidewire introducer is used, remove according to manufacturer's instructions.
- Slide torque device over proximal end of guidewire per Figure 1, and secure.

Steam Shaping

Excelsior® XT-17™ Microcatheter is packaged with a steam shaping mandrel for the physician to add shape proximal to the distal tip if desired.

To maintain the integrity of the microcatheter inner diameter, these steam shaping instructions should be strictly followed.

Warning: The shaping mandrel is not intended for use inside the human body.

- Remove shaping mandrel from packaging and insert mandrel's entire working length into microcatheter's distal lumen. Flush or dip microcatheter tip into heparinized saline to aid mandrel movement.
- Bend working length of shaping mandrel, proximal to the distal tip, to desired shape.
- Shape microcatheter by holding mandrel/microcatheter assembly no closer than 2.54 cm (1 in) from steam source for approximately 10 seconds.

Caution: Do not position microcatheter closer than 2.54 cm (1 in) from the steam source. Damage to the microcatheter may result.

- Remove mandrel from microcatheter and discard. Multiple shaping is not recommended.

Directions For Use

- Place appropriate guide catheter according to the manufacturer's instructions.
- Open RHV thumb-screw and carefully introduce microcatheter and guidewire assembly through the guide catheter RHV using the peel away introducer. Tighten the O-ring valve around the microcatheter just enough to prevent back flow, but not so tightly as to inhibit microcatheter advancement. Remove the peel away introducer from the microcatheter shaft by retracting the introducer from the RHV and pulling the tab.

Caution: Excessive tightening of a hemostatic valve onto the microcatheter shaft may result in damage to the microcatheter. Removing the peel away introducer without a guidewire inserted in the microcatheter lumen might result in damage to the microcatheter shaft.

- Advance the microcatheter until the fluoro saver marker reaches the guide catheter RHV.
- Under fluoroscopic guidance advance microcatheter tip by carefully rotating torque device on proximal end of guidewire. Hold proximal end of guidewire and microcatheter in a straight line to enhance transmission of rotational energy from proximal to distal end.
- Advance the guidewire a short distance, then alternately track the microcatheter over the guidewire until selected site is reached. The guidewire may also be slowly retracted during microcatheter advancement. This step is best accomplished by two people.

Caution: To facilitate microcatheter handling, the proximal portion of the microcatheter does not have the hydrophilic surface. Greater resistance may be encountered when this section of the microcatheter is advanced into the RHV.

- To infuse, completely remove the guidewire from the microcatheter, connect 3 cm³ (3 cc) syringes, and infuse as required. Monitor infusion pressures with a syringe manometer device whenever practical.

Warning: Discontinue use of microcatheter for infusion if increased resistance is noted. Resistance indicates possible blockage. Remove and replace blocked microcatheter immediately. DO NOT attempt to clear blockage by over-pressurization. Doing so may cause the microcatheter to rupture, resulting in vascular damage or patient injury.

Table 2. Flow Rate Table

Microcatheter	Shaft Length (cm)	Microcatheter ID (mm/in)	Approximate Dead Space (cm³)	Approximate Flow Rates at 2,070 kPa (300 psi) cm³/sec		
				Water (cm³/sec)	60% Contrast (cm³/sec)	76% Contrast (cm³/sec)
Excelsior XT-17	150	0.43 mm (0.017 in)	0.30	0.83	0.20	0.09

Warning: Do not exceed 2,070 kPa (300 psi) infusion pressure. Excessive pressure could dislodge a clot, causing thromboemboli, or could result in a ruptured microcatheter or severed tip, causing vessel injury.

WARRANTY

Stryker Neurovascular warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Stryker Neurovascular's control directly affect the instrument and the results obtained from its use. Stryker Neurovascular's obligation under this warranty is limited to the repair or replacement of this instrument and Stryker Neurovascular shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Stryker Neurovascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Stryker Neurovascular assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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Copyright © 2015 Stryker
Stryker Australia Pty Ltd