Surpass Evolve[™]

Flow Diverter System

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	13
Gebrauchsanweisung	19
Istruzioni per l'uso	25
Gebruiksaanwijzing	31
Instruções de Utilização	37

Surpass Evolve™

Flow Diverter System

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Stryker Neurovascular representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Surpass Evolve Flow Diverter System is comprised of a self-expandable braided implant preloaded on a delivery wire, housed inside an introducer sheath. The system consists of the following components:

- Surpass Evolve Flow Diverter (implant)
- Delivery wire and introducer sheath
- Optional torque device

Surpass Evolve Flow Diverter

The braided Surpass Evolve Flow Diverter is the implant portion of the system. It is shipped pre-loaded on the delivery wire, inside the introducer sheath. These devices come in different diameter and length combinations ranging from 2.5mm to 5.0mm in diameter and from 12mm to 40mm in length. Interwoven with the Surpass Evolve Flow Diverter cobalt chromium wires are platinum-tungsten wires for visualization under fluoroscopy.


Once released from the constraint of the microcatheter into the vessel, the implant expands to the vessel lumen diameter and can be recaptured and repositioned prior to full deployment, if needed. In its expanded shape, the implant diverts the blood flow away from the aneurysm.

Delivery Wire and Introducer Sheath

The implant is pre-loaded on the delivery wire and protected by an introducer sheath. The flow diverter is located between the distal and proximal markers on the delivery wire, which signal the distal and proximal ends of the implant. The petal is designed to protect the distal end of the implant as it is advanced through the microcatheter. The resheath pad allows the user to resheath the flow diverter back into the microcatheter. The resheath marker provides the operator fluoroscopic visualization of the resheathing limit of the implant. The distal 2cm of the delivery wire, the distal tip, is radiopaque.

The introducer sheath consists of a polymer tube with a tapered distal end. It is designed to protect the implant until it is transferred into the microcatheter.

Figure 1. Surpass Evolve Flow Diverter System

- | | | |
|--------------------|--------------------|----------------------|
| 1. Delivery Wire | 4. Resheath Marker | 7. Distal Marker |
| 2. Proximal Marker | 5. Flow Diverter | 8. Distal Tip |
| 3. Resheath Pad | 6. Petal | 9. Introducer Sheath |

Optional Torque Device

The physician may attach the torque device to the proximal end of the delivery wire, which may facilitate advancement through the microcatheter and deployment. However, the delivery wire is not designed to be torqued.

Accessories (Not Provided)

Standard interventional devices, including:

- Appropriately sized Sheath or Guide Catheter
- Appropriately sized Intermediate Catheter (recommended)
- Stryker Neurovascular Excelsior® XT-27® Standard Straight (0.027 in ID) Microcatheter with a 135cm or 150cm length
- Appropriately sized Guidewire
- Rotating hemostatic valves \geq 4F (9F recommended)
- Sterile heparinized saline solution

The Surpass Evolve Flow Diverter System is designed specifically for use with a Stryker Neurovascular Excelsior XT-27 Standard Straight Microcatheter.

Compatibility with other microcatheters has not been established.

INTENDED USE/INDICATIONS FOR USE

The Surpass Evolve Flow Diverter System is indicated for use for the treatment of saccular or fusiform intracranial aneurysms arising from a parent vessel with a diameter \geq 2.0mm and \leq 5.0mm.

CONTRAINDICATIONS

The Surpass Evolve Flow Diverter is contraindicated in the following patient types:

- Patients in whom the parent vessel size does not fall within the indicated range.
- Patients in whom antiplatelet and/or anticoagulation therapy (e.g., aspirin and clopidogrel) is contraindicated.

- Patients who have not received dual anti-platelet agents prior to the procedure.
- Patients with an active bacterial infection.
- Patients in whom the angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as:
 - Severe intracranial vessel tortuosity or stenosis; and/or
 - Intracranial vasospasm not responsive to medical therapy.

WARNINGS

- This device should only be used with adequate fluoroscopic guidance by physicians who have received appropriate training in interventional neuroradiology or interventional radiology, and training on the use of this device as established by Stryker Neurovascular.
 - Persons allergic to nickel, cobalt chromium or platinum tungsten metal may suffer an allergic response to this system.
 - Appropriate anti-platelet and anti-coagulation therapy should be employed in accordance with standard medical practice.
 - A thrombosing aneurysm may aggravate pre-existing or cause new symptoms of mass effect and may require medical therapy.
 - Replace the system if any defects are observed during use.
 - Purge the entire system carefully to avoid the accidental introduction of air into the patient.
 - At any time during the procedure discontinue use of the system if excessive resistance is encountered. Movement of the system against resistance may result in damage to the vessel, a system component, or the patient.
 - Do not resheath the implant into the introducer sheath once transferred into the microcatheter as it may result in damage to the implant wires.
 - It is important for the Surpass Evolve™ Flow Diverter System to be used with the Stryker Neurovascular Excelsior® XT-27® Standard Straight Microcatheter with a 135cm length (Model UPN M003XT2735810) or 150cm length (Model UPN M003XT2750810).
 - Repositioning of the Surpass Evolve Flow Diverter in the parent vessel without fully retrieving the device is not advised since it could cause vessel damage and/or perforation. After full deployment, do not re-position.
 - [Clinical Warning] Do not use Surpass Evolve Flow Diverter for ruptured intracranial aneurysms.
 - [Clinical Warning] A decrease in effectiveness has been observed in subjects aged > 65 years old, subjects with history of smoking and history of prior non-target intracranial aneurysm treated.
 - [Clinical Warning] Judicious patient selection is important. Patients who fall outside the therapeutic range for antiplatelet testing or at the lower limits of platelet inhibition have an increased risk of developing stent thrombosis, even with additional doses of antiplatelet medication. Another effective anti-platelet agent should be considered.
 - [Clinical Warning] Placement of multiple Surpass Evolve Flow Diverters may increase the risk of ischemic complications.
 - [Clinical Warning] Delayed aneurysm rupture may occur with large and giant intracranial aneurysms.
-

PRECAUTIONS

- This device has not been evaluated for pediatric use.
- The safety and effectiveness of the device has not been established in the treatment of small and medium wide-neck intracranial aneurysms.
- Experience with endovascular implants indicates that there is a risk of stenosis. Subsequent stenosis may require dilatation of the vessel segment containing the device. The risks and long-term outcome following dilatation of endothelialized devices is unknown at present.
- Select a device length that is at least 10mm longer than the intracranial aneurysm neck to maintain a minimum of 5mm on either side of intracranial aneurysm neck. Confirm the device labeling reflects the desired size of the target vessel where the device is to be used.
- Use care when removing the device from its hoopcard so as not to eject the implant from the introducer sheath.
- Do not expose the system to organic solvents (e.g., alcohol) as it may result in damage to the system.
- Do not torque or rotate the delivery wire as it may result in system damage.
- Do not attempt to partially deploy and recapture the implant more than three times as it may result in a loss of system performance.
- Use care when crossing the deployed implant with guidewires or other accessory devices as it may cause the implant to migrate distally.
- Operators should take all necessary precautions to limit X-radiation dose to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.
- Lower intracranial aneurysm occlusion rates may be associated with giant intracranial aneurysms (>25mm).
- Lower intracranial aneurysm occlusion rates may be associated with implants that are not fully apposed to the vessel wall.
- Do not deploy Surpass Evolve Flow Diverters in parallel (side by side).

ADVERSE EVENTS

Risks that may be associated with the use of the Surpass Evolve Flow Diverter System in the intracranial arteries include:

- Adverse reaction to anesthesia, contrast or antiplatelet/anticoagulation agents
- Allergic reaction
- Aphasia
- Cardiac Arrhythmia
- Confusion, coma, change in mental status
- Cranial neuropathy/neurological deficits
- Death
- Device migration, fracture, misplacement
- Dissection of the parent artery
- Embolism (air, clots, device fragments)
- Groin injury (bleeding, pain, vessel/nerve damage)
- Headache
- Hemiplegia
- Hydrocephalus
- Implant or parent vessel stenosis
- Implant thrombosis/occlusion
- Infection
- Intracerebral bleeding

- Mass effect
- Myocardial infarction
- Neurological deficits
- Perforation or rupture of aneurysm
- Perforation of parent artery
- Progressive neurologic symptoms related to intracranial aneurysm (IA)
- Pseudoaneurysm formation
- Reaction to radiation exposure (i.e. alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, delayed neoplasia)
- Renal failure
- Retroperitoneal hematoma
- Seizure
- Stroke
- Subarachnoid hemorrhage
- Thromboembolism from device
- Thrombosis of parent artery or branch vessel
- Transient ischemic attack (TIA)
- Vasospasm

Risks that are eye related with the use of the Surpass Evolve™ Flow Diverter may include:

- Amaurosis fugax/transient blindness
- Blindness
- Diplopia
- Reduced visual acuity/field
- Retinal artery occlusion
- Retinal ischemia
- Retinal infarction
- Vision impairment

Please notify your Stryker Neurovascular representative immediately if a device malfunctions or patient complication or injury is experienced or suspected. Please make every attempt to retain any suspect device, its associated components and their packaging for return to Stryker Neurovascular.

MRI SAFETY INFORMATION

Non-clinical testing demonstrated that the Surpass Evolve Flow Diverter is MR Conditional for single and overlapping up to 60mm in length.

A patient with the Surpass Evolve Flow Diverter can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only
- Maximum spatial gradient magnetic field of 3,000 Gauss/cm (30 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg in the Normal Operating Mode

Under the scan conditions defined, the Surpass Evolve Flow Diverter is expected to produce a maximum temperature rise of 3.0°C after 15 minutes of continuous scanning (i.e., per pulse sequence). In non-clinical testing, the image artifact caused by the Surpass Evolve Flow Diverter extends approximately 10mm from this implant when imaged using a gradient echo pulse sequence and a 3 Tesla MR system. The lumen cannot be visualized on gradient echo or T1-weighted, spin echo pulse sequences.

Caution: The Surpass Evolve Flow Diverter may create local field inhomogeneity and susceptibility artifacts during magnetic resonance angiography (MRA), which may degrade the diagnostic quality to assess effective intracranial aneurysm treatment.

The health state of the patient or the presence of other implants may require reduction of the MRI limits.

HOW SUPPLIED

Stryker Neurovascular products are single use only, sterile and non-pyrogenic in unopened packaging that is designed to maintain sterility unless the primary product pouch has been opened or damaged.

WARNINGS

- Do not use if sterile barrier or package is opened or damaged. If damage is found, call your Stryker Neurovascular representative.
- Do not use if labeling is incomplete or illegible.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

Handling and Storage

Store in a cool, dry, dark place.

PREPARATIONS FOR USE

Surpass Evolve Flow Diverter System Selection

Appropriate selection of the device is important for patient safety. In order to choose the optimal device for any given intracranial aneurysm neck length and parent artery diameter, examine preprocedure angiograms.

DIRECTIONS FOR USE

Flow Diverter Selection and Preparation

1. Select an appropriate Surpass Evolve Flow Diverter System based on the sizing recommendations in **Table 1**.

Table 1. Device Sizes

Device Diameter (mm)	Device Length (mm)	Unconstrained Implant Outer Diameter (mm)	Recommended minimum vessel diameter (mm)	Recommended maximum vessel diameter (mm)
2.5	12, 13, 14, 15, 16, 17, 18, 19, 20	>2.5	≥2.0	≤2.5
3.25	12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 25, 30	>3.25	>2.5	≤3.25
4.0	12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 25, 30, 40	>4.0	>3.25	≤4.0
4.5	12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 25, 30, 40	>4.5	>4.0	≤4.5
5.0	12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 25, 30, 40	>5.0	>4.5	≤5.0

Note: The device diameter should be chosen based on the measurements of the vessel size proximal to the intracranial aneurysm neck.

2. Choose a Surpass Evolve™ Flow Diverter System with a labeled length that is at least 10mm longer than the intracranial aneurysm neck (to maintain a minimum of 5mm on either side of the neck) and allows for the proximal and distal ends of the implant to land in a straight section of the vessel.

Note: The implant changes in length during deployment, referred to as foreshortening. Take this into account when deploying (Table 2). The implant deployed length will never be shorter than its labeled length, in its indicated vessel diameter, even after foreshortening.

Table 2. Device Deployed Length Changes*

Device Diameter (mm)	Average decrease in length from inside Delivery System to fully expanded (%)
2.5	38
3.25	35
4.0	49
4.5	52
5.0	58

*Based on mathematical calculations

Flow Diverter Positioning

1. Using conventional catheterization technique, place the microcatheter tip at least 20mm past the distal neck of the intracranial aneurysm. It is recommended to use a heparinized saline drip to continuously flush the microcatheter during the procedure.
2. Remove the guidewire from the microcatheter.
3. Remove the Surpass Evolve Flow Diverter System packaging hoop from the pouch and pull the proximal end of the introducer sheath from the white clip on the packaging hoop. Carefully remove system from the packaging hoop.

Note: Hold the introducer sheath and delivery wire together when removing the system from the packaging hoop to avoid premature deployment.

4. Inspect the system for any damage. Confirm that the implant and delivery wire are properly pre-loaded in the distal tip of the introducer sheath.
5. Partially insert the distal tip of the introducer sheath into the microcatheter RHV and tighten the RHV.

Note: Under-tightening or over-tightening the RHV may result in inadequate flushing.

6. Using an appropriately pressurized sterile flush or syringe, confirm that fluid exits the proximal end of the introducer sheath prior to advancing the system into the microcatheter.

Warning: Confirm that there are no air bubbles in the system.

7. Loosen the RHV, advance the introducer sheath into the RHV and visually confirm the tip of the sheath is seated deeply in the hub of the microcatheter. Tighten the RHV firmly to secure the introducer sheath to the hub of the microcatheter.

8. Advance the delivery wire to transfer the implant from the introducer sheath into the microcatheter.
9. Continue to advance the system until the distal edge of the fluoro saver marker enters the proximal end of the introducer sheath.
10. Loosen the RHV on the microcatheter and remove the introducer sheath, holding the delivery wire in place. Re-tighten the RHV as needed. If desired, place a torque device on the proximal end of delivery wire to facilitate advancement through the microcatheter and deployment.

Caution: Use the torque device when advancing the delivery wire only. Do not use the torque device to torque the delivery wire as it is not designed to be torqued and it may result in system damage.

11. Slowly advance the delivery wire until the distal edge of the fluoro saver marker reaches the microcatheter RHV. The fluoro saver marker is located 125cm from the distal tip of the delivery wire.

Warning: The fluoro saver marker is only compatible with microcatheters with a minimum length of 135cm. Usage of shorter microcatheters could result in inadequate working length for the operator and unintentional deployment of the implant.

12. Under fluoroscopic visualization, continue to advance the system through the microcatheter until the distal radiopaque marker on the delivery wire is at least 5mm beyond the neck of the intracranial aneurysm.
13. Confirm the position of the implant by visualizing the distal and proximal radiopaque markers on the delivery wire. Ensure the distal end of the implant is at least 5mm beyond the neck of the intracranial aneurysm.

Note: Pulling back on the microcatheter to make final adjustments will ensure that slack has been removed from the system prior to deployment.

14. Adjust RHV tightness just enough to allow free movement between the microcatheter and delivery wire. The Surpass Evolve Flow Diverter is now ready to be deployed.

Flow Diverter Deployment

1. Begin deployment by retracting the microcatheter while putting slight forward pressure on the delivery wire until the implant begins to exit the microcatheter and starts apposing to the vessel wall.

Warning: Once the implant has begun to appose the vessel walls, do not attempt to move the implant more distally as damage to the implant ends may occur.

2. A combination of forward pressure (pushing) on the delivery wire and/or retracting the microcatheter (unsheathing) should enable the operator to keep the distal end of the implant aligned at the desired landing zone as well as achieve optimal vessel wall apposition.

Note: During deployment, wall apposition may be improved by:

- Resheathing the implant and/or manipulation of the microcatheter by locking down the delivery wire using the RHV and gently advancing both the delivery wire and the microcatheter as a system (loading); and/or
- Maintaining the microcatheter tip within the center of the parent vessel.

3. Resheathing instructions: If positioning is not satisfactory, the implant can be recaptured and repositioned as long as the distal marker on the microcatheter is not proximal to the resheathing marker on the delivery wire. To recapture the implant, simultaneously advance the microcatheter while gently pulling back on the delivery wire.

Caution: Do not attempt to partially deploy and recapture the implant more than three times or a loss of resheath performance may occur.

Caution: Do not apply excessive force if significant resistance is experienced while attempting to recapture the implant or a failure of the device may occur.

4. After about 10mm of the implant has been deployed, its distal end will begin to flare and appose to the vessel walls.
5. After the distal end of the implant has successfully expanded and begun apposing to the vessel walls, deploy the remainder of the implant by alternately pushing the delivery wire and retracting the microcatheter. Carefully monitor the delivery wire distal tip under fluoroscopy during deployment of the implant. Once the operator has passed the resheath marker, the implant should be unsheathed without pushing on the delivery wire.
6. After the entire implant has deployed, confirm full expansion under fluoroscopy and ensure that it has completely apposed the vessel wall. If the device is not fully apposed, consider using a balloon catheter to fully open it.
7. Advance the microcatheter over the delivery wire through the deployed implant, being careful not to dislodge the implant, until the radiopaque marker on the microcatheter tip is distal to the deployed implant. Retract the delivery wire into the microcatheter.
8. Verify that the device has remained patent and properly positioned.
9. After completing the procedure, dispose of all devices in accordance with hospital, administrative and/or local government policy.

Frequently Asked Questions (FAQ)

- Q. If excessive friction is experienced during the insertion of the delivery system into the microcatheter, what should I do?
 - A. Carefully remove the entire system simultaneously (microcatheter and delivery system).
- Q. Can I retrieve the implant if the distal end has deployed at an undesirable location?
 - A. Yes. A partially deployed implant can be resheathed. Carefully pull back the delivery wire while simultaneously advancing the microcatheter until the implant is secured inside the microcatheter. Then, if there is no resistance, simultaneously remove the entire system (microcatheter and delivery system).
- Q. Can I retrieve a fully deployed implant?
 - A. Once fully deployed, the implant cannot be removed. A second implant can be deployed, if needed.
- Q. Can I place a second implant inside another implant?
 - A. Yes. A second implant can be placed inside another. After placing the first implant, advance the microcatheter over the delivery wire while keeping the delivery wire across the neck of the intracranial aneurysm. Position the microcatheter at the desired location and retrieve the delivery wire. Select a new appropriate implant and deploy it as normal.

Warning: Placement of multiple Surpass Evolve™ Flow Diverters may increase the risk of ischemic complications.

- Q. If there is a difference between the proximal and distal vessel diameters, which Surpass Evolve Flow Diverter System diameter do I choose?
 - A. Choose a system that matches the larger (typically proximal) vessel diameter to ensure proper anchoring.

WARRANTY

Stryker Neurovascular warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Stryker Neurovascular's control directly affect the instrument and the results obtained from its use. Stryker Neurovascular's obligation under this warranty is limited to the repair or replacement of this instrument and Stryker Neurovascular shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Stryker Neurovascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Stryker Neurovascular assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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Copyright © 2018 Stryker
Stryker Australia Pty Ltd